ENGLISH

Use the Woller as a Heating Pad:

- Relaxes and reduces pain due to abdominal cramps and congestion.
- For providing relaxing heat and pressure stimulation when the baby cries and is irritable.
- Facilitates breastfeeding. Place a warm Woller on the breast.
- Pain relief for engorgement. Place a cold Woller on the breast.
- In case of ligament pain.
- For lower back pain during pregnancy (LBP)
- For preheating the bed, towel and pram cradle.

User Instructions KipKep Woller:

- To be used as a heating pad and cold compress.
- The shape of the Woller makes it ideal for the baby's belly in case of colic. Lay your baby on his/her back and place the Woller softly on the belly. The weight is distributed along the sides of the belly.
- For safety reasons, always read and follow these instructions before usage. Do not discard them.
 - The Woller should always feel lukewarm to the touch, never hot!
 - If the microwave has a grill, always remove it!
 - Do not exceed the recommended heating time and only reheat the pad when it has cooled down to room temperature.
 - Always heat the Woller on the rotating glass plate, preferably on a clean microwave-proof (ceramic or glass) plate. Always consult your microwave's manual.
 - Ensure that your microwave is free of food scraps. Food scraps can become very hot when reheated and might scorch patches of the Woller's fabric.
 - Allow the Woller to cool for 1 minute before applying it to a body. After this minute, always check the temperature on the inside of your wrist or in the cavity of your arm if the warmed up pillow is suitable for babies (up to 40°C).
 - Always heat the pad without a cover.
 - Never use the heating pad directly on bare skin.
 - Wash the cover before using it for the first time. In order to clean it, you can wash the outer cover at 40°C degrees.
 - Place the Woller on the belly as shown in the illustration 4.

Heating:

- (Ill. 1) Remove the pillow out of the cover.
- Place the pad gently in the microwave, ensuring that nothing interferes with the movement of the turntable. If your microwave does not have a turntable, you should place the pad on a microwave-safe plate and turn the pad halfway through the heating time to ensure that the pad is thoroughly heated. Take one moment somewhere in between to shake the Woller back and forth.
- If the microwave has a grill, always remove it before heating the pad and make sure the pad is placed on a clean surface. Never lay it on a metal surface and make sure it cannot touch the walls during heating. If preferred, fold the pillow in half; seeds evenly distributed. If your microwave doesn't have a grill, you should place the pad on a microwave-proof plate/bowl and turn it halfway through heating to ensure that the product is heated evenly.
- (Ill. 1) Heat the pad without the cover for 30 seconds in a microwave at 600 Watts. Repeat with 10 sec. until the desired temperature is reached. Do not heat the pad for longer than 60 seconds.
- Carefully remove the pad from the microwave to prevent burns. Always check the temperature as described above. Comfortably warm is warm enough for a baby.
- The pad should feel lukewarm, not hot, and never too hot to touch. Overheating can cause a fire or scorching. Overheating can damage the pad as well as the microwave.
- The microwave heating process slowly distributes the heat from the inside to the outside. Therefore, it will take some time before the pad feels warm after removing it from the microwave.
- Put the cushion back in the cover.
- (Ill. 2) Check the temperature of the Woller 1 minute after heating it to verify whether it has cooled down sufficiently to be used.
- (Ill. 3) Always check the temperature as described above.
- Place the pad back in the cover.
- Do not leave the Woller in the microwave for longer than the recommended time to prevent overheating! Do not leave the pad unattended during heating.
- (Ill. 4) Place the Woller in the cover on the baby's covered tummy. Always use under supervision.
- The pad can also be used as a cold compress by storing it in the refrigerator or freezer in a clean freezer bag.

Cooling:

- For use as a cold compress.
- Cool the cushion with the cover in a clean freezer bag in a refrigerator or freezer.
 - Keep the Woller in the refrigerator/freezer for at least 30 minutes.
 - Always use a clean freezer bag when cooling the Woller to prevent contact with food and/or ice.

How to store:

- Store in a cool, dry place.
- Do not store in direct sunlight.

Warning:

- Do not use the pad in a conventional or hot air oven.
- Do not allow any contact between the Woller and metal when it is being heated, for example the grate, the wall, or the bottom of the microwave. This damages the Woller.
- Always check the temperature of the pad before usage.
- Do not use the pad on sensitive skin.
- Do not use if symptoms of skin irritation occur.
- Do not use in case of symptoms of inflammation.
- Do not allow the pad to become wet. If damp, let the heating pad dry in a warm spot. If the pad has become

- wet, it can no longer be used.
- Do not use while sleeping or on others that are sleeping.
- The heating pad itself CANNOT be washed. However, the cover is washable.
- The pad is no longer useable after overheating (either the seeds char and may give off an unpleasant smell or the pad or cover is damaged by burns)
- Only use the pad on the baby's lower body, and never use it without supervision.
- KipKep is not liable for any damage whatsoever resulting from incorrect usage or failure to comply with the user/manufacturer instructions of the pad.
- Any serious incident that occurs during use of the Woller should be reported to the manufacturer (KipKep) and the competent authority of the country in which the user/patient resides.
- Always follow your microwave manufacturer's instructions and check the wattage per setting of your microwave.
- Composition of the Woller. Filling: linseed | Pillow: 100% cotton | Cover: varies, see label in cover. Always use the Woller with the original KipKep cover: Woller Cover.

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KipKep B.V.